


Chemistry 1105 Lab: Fatty Acids

Goals:

- 1. Introduction to Titration.**
- 2. Prepare and standardize 0.1 M NaOH.**
- 3. Determine the Percentage of Acetic Acid(weight/volume) in Vinegar.**

Acid-Base Titration:

A titration is a process in which a controlled volume of one reagent(titrant) is added to a known amount or volume of a second reagent until a complete reaction is observed.



Preparation of 0.1 M NaOH:

Preparation of 150. mL of a solution of 0.1 M NaOH by dilution of a 1 M NaOH stock solution.

$$M_1 \times V_1 = M_2 \times V_2$$


M_1 : Molarity of stock solution

V_1 : Volume of stock solution

M_2 : Molarity of dilute solution

V_2 : Volume of dilute solution

Standardization of 0.1 M NaOH:



? M
mL

g
↓

moles $\text{KHC}_8\text{H}_4\text{O}_4$

moles $\text{KHC}_8\text{H}_4\text{O}_4 = \text{moles NaOH}$ at
equivalence point

Determination of the Endpoint/Equivalence

Point:

Equivalence point determined using acid-base indicator.



Equivalence point is volume of base where the moles base = moles acid.


Endpoint. Volume of base that turns indicator color.

Determining Percent Acetic Acid in Vinegar:

$$\text{Percent Acetic Acid}\left(\frac{w}{v}\right) = \frac{\text{mass of acetic acid(g)}}{100.00 \text{ mL of solution}}$$

**NOTE: Can not analyse unaltered vinegar.
Will perform 1/10th dilution.**

Will analyse 25.00 mL of this diluted vinegar.



(now known)

0.1 M
mL



moles NaOH

=

M = moles/volume(L)
25.00 mL DILUTE



moles CH₃COOH

Finding the Percent Acetic Acid:

Percent Acetic acid is the mass of acetic acid(g) in 100.00 mL.

$$\text{Percent Acetic Acid} \left(\frac{\text{w}}{\text{v}} \right) = \text{Molarity} \left(\frac{\text{moles}}{\text{L}} \right) \times \frac{60.05 \text{ g acetic acid}}{\text{mole}} \times 0.10000 \text{ L}$$

Expected value: 7.0% w/v