

Chemistry 1104 Lab: Volumetric

Goals:

- 1. Introduction to Volumetric Glassware.**
- 2. Precision of Volumetric Glassware.**
- 3. Introduction to Concentration
Molarity**
- 4. Preparation of a Solution of known
Molarity from a) solid b) dilution.**

Volumetric Glassware: Graduated Cylinder

TC: To Contain
Least Precise

Sizes

10.0 mL


25.0 mL

50.0 mL

100.0 mL


Precision

± 0.1 mL




Reading a Graduated Cylinder:

TC: To Contain



34.2 mL



34.2 mL

Volumetric Flask

TC: to contain

Used to prepare solutions of known concentration.

Sizes

25.00 mL

50.00 mL

100.00 mL

250.00 mL

Precision

± 0.01 mL



Volumetric Pipet

TD: To Deliver

Used to dispense specific volumes of liquid.

**Sizes: 1.00 mL, 2.00 mL,
5.00 mL, 10.00 mL, 20.00 mL,
25.00 mL**

**Precision
 ± 0.01 mL**




Buret

TD: To Deliver

Used to dispense known amounts of liquid of various volume.

Size: 50 mL max capacity typically. Precision ± 0.01 mL

**When Full Volume recorded
0.00 mL**




Reading a Buret:

When full the volume recorded with a buret is **0.00 mL**.

Buret TD.

Full buret not recorded as 50 mL. WRONG!



8.20 mL

Part A:

**Determining volume of a test tube using a
Graduated cylinder.**

Part B:

**Determining volume of same test tube using a
buret.**

Part C:

Determining the precision of a volumetric pipet by measuring the mass of water Dispensed by the pipet.

Record mass of small beaker before and after. Recommend using 10 mL pipet. Volume recorded as 10.00 mL.

NOTE: mass(g) of water = volume(mL)

Part D:


Introduction to Molarity.

$$\text{Molarity(M)} = \frac{\text{moles of solute}}{\text{volume of solution(L)}}$$

**For example a 4.00 M HCl solution indicates
That there is 4.00 mol HCl for every 1 L of
solution.**

Part D: cont...

Preparing a solution using a volumetric flask and a known mass of a solid(solute).



MW = 249.68 g/mol

Assigned a Molarity.

Assigned Volume: 100.00 mL



Part D: cont...

The mass of solute depends on Molarity and the volume of solution needed.

$$\mathbf{100.00\text{ mL} = 0.10000\text{ L}}$$

$$0.10000\text{ L} \times \frac{\text{moles of solute}}{\text{volume of solution(L)}} (\text{M}) = \text{moles solute}$$

$$\text{moles solute} \times \frac{\text{mass of solute(g)}}{\text{mole solute}} (\text{Molar Mass}) = \text{g of solute}$$

Part D: cont...

Weigh out desired mass of solute.

Record mass used.

Mass of solute used → moles solute

$$\text{Molarity(M)} = \frac{\text{moles of solute}}{0.10000(\text{L})}$$

Part D: cont...

Use spec20 to determine if solution prepared correctly. Spec20 measures the amount of light absorbed by the aqueous solution.

Record Absorbance.

Part E:

Preparation of a solution of known concentration by dilution of a stock solution of known concentration.

$$M_1 \times V_1 = M_2 \times V_2$$

M_1 : Molarity of stock solution

V_1 : Volume of stock solution

M_2 : Molarity of dilute solution

V_2 : Volume of dilute solution