## **Chemistry 1105 Lab: ASA Titration**

#### **Goals:**

- **1. Technique of Titration.**
- 2. Prepare and standardize 0.1 M NaOH.
- **3.** Determine %ASA in the sample you synthesized in previous experiment.

#### **Acid-Base Titration:**

A titration is a process in which a controlled volume of one reagent(titrant) is added to a known amount or volume of a second reagent until a complete reaction is observed.

**Base** + Acid  $\rightarrow$  Salt + Water

## Preparation of 0.1 M NaOH: Preparation of 200. mL of a solution of 0.1 M NaOH by dilution of a 1 M NaOH stock solution.

 $\mathbf{M}_1 \!\times\! \mathbf{V}_1 \!=\! \mathbf{M}_2 \!\times\! \mathbf{V}_2$ 

**M<sub>1</sub>: Molarity of stock solution** 

- **V<sub>1</sub>: Volume of stock solution**
- M<sub>2</sub>: Molarity of dilute solution
- **V<sub>2</sub>: Volume of dilute solution**

#### **Standardization of 0.1 M NaOH:**

## Base + Acid $\rightarrow$ Salt + Water NaOH + KHC<sub>8</sub>H<sub>4</sub>O<sub>4</sub> $\rightarrow$ KNaC<sub>8</sub>H<sub>4</sub>O + H<sub>2</sub>O ? M g mL $\downarrow$ moles KHC<sub>8</sub>H<sub>4</sub>O<sub>4</sub>

#### moles KHC<sub>8</sub>H<sub>4</sub>O<sub>4</sub> = moles NaOH at equivalence point

## **Determination of the Endpoint/Equivalence Point:**


## **Equivalence point determined using acidbase indicator.**



Equivalence point is volume of base where the moles base = moles acid.

**Endpoint. Volume of base that turns indicator color.** 

#### **Determining the %ASA of Sample:**



# o-C<sub>6</sub>H<sub>4</sub>(OCOCH<sub>3</sub>)COOH

#### C<sub>6</sub>H<sub>4</sub> aromatic ring o-ortho (position of substituents)

$$\% ASA = \frac{mass of ASA(g) from titration}{mass of sample(g) weighed} \times 100\%$$

## Base + Acid $\rightarrow$ Salt + Water NaOH + 0-C<sub>6</sub>H<sub>4</sub>(OCOCH<sub>3</sub>)COOH $\rightarrow$ 0-C<sub>6</sub>H<sub>4</sub>(OCOCH<sub>3</sub>)COONa + H<sub>2</sub>O

