



# Chemistry 2201 Exp: GASES

**Determining the Heat Capacity of a gas at constant pressure( $C_p$ ) and constant volume ( $C_v$ ).**

**Acieved by measuring 3 pressures of a fixed amount of gas.**

# Clément and Desormes Method





**NOTE: Gas is compressed.**

**$P_1$ ,  $T_1$ ,  $V_1$ : Initial pressure, temperature and volume of the gas.**

**$P_1$  greater than atmospheric pressure( $P_2$ ).**

**$V_1$  is the compressed volume of gas under study. Not volume of container.**

**Step 1:**

**Gas is allowed to expand from  $P_1$  to atmospheric pressure( $P_2$ ).**

**As gas expands it cools from  $T_1$ (room temperature) to  $T_2$ . Occupies volume of vessel( $V_2$ ).**

**Heat flows into vessel. Volume constant. Pressure increases to  $P_3$  as gas heats back to room temperature.**

**Assume gas behaves Ideally.**

**For an ideal gas the internal energy only depends on temperature.**

**By measuring three pressures  $P_1$ ,  $P_2$ , and  $P_3$  the heat capacity ratio( $\gamma$ ) is determined.**

$$\gamma = \frac{\ln\left(\frac{P_1}{P_2}\right)}{\ln\left(\frac{P_1}{P_3}\right)} = \frac{\ln P_1 - \ln P_2}{\ln P_1 - \ln P_3}$$

**Can determine  $C_P$  and  $C_V$  from  $\gamma$ .**

$$\gamma = \frac{C_P}{C_V} = \frac{C_V + R}{C_V} \quad C_P = C_V + R$$

**Pressure measured using dibutyl phthalate manometer. Measures pressure difference in vessel to atmospheric pressure.**

$$\text{cm DP} \times \frac{\text{density DP}}{\text{density Hg}} = \text{cm Hg}$$

**Density: dibutyl phthalate =  $1.046 \text{ g}\cdot\text{cm}^{-3}$   
mercury =  $13.55 \text{ g}\cdot\text{cm}^{-3}$**